

# Whatman™ 934-AH™ Glass Microfiber filters

Whatman 934-AH Glass Microfiber filters (Fig 1) are suitable for suspended solids analysis. These filters also permit rapid passage of very large sample volumes without blockage for the collection of suspended solids in potable, natural, and industrial waters. Specified in standard methods\* for determining total suspended solids in water, Whatman 934-AH Glass Microfiber Filters are excellent for the removal of turbidity and filtration of bacterial cultures. The fine particle retention of this popular grade is superior to any cellulose filter of similar rating and is capable of high retention efficiency at high flow rates.

934-AH is now available in a RTU (ready to use) format for your Total Suspended Solids.

## Features and benefits

- **Excellent fine particle retention:** allows for high particle retention efficiency when filtering large volumes.
- **Binder-free, high-grade borosilicate glass microfiber:** enables use with temperatures up to 550°C.
- **Available in a wide range of sizes:** in circle and sheet formats.
- **High loading capacity:** capable of processing very turbid streams.
- **Low fiber shedding:** improves quality assurance of test results.

## Applications

- Water pollution monitoring
- Air pollution monitoring
- Cell harvesting
- Liquid scintillation counting

\*The most widely used filter for establishing water quality in suspended solids content and related measure (Standard Method 2540D and EPA Method 160.2). Total suspended (non-filterable) solids are defined as those which are retained by a standard glass microfiber filter without organic binder.



Fig 1. Whatman 934-AH Glass Microfiber filters are designed for the analysis of suspended solids.




Fig 2. A generalized protocol for the analysis of suspended solids with Whatman 934-AH Glass Microfiber filters.

## Typical data

Weight (g/m <sup>2</sup> )	64
Thickness (µm)	435
Particle retention (µm)	1.5
Filtration speed (Sec/100 ml)*	47
DOP penetration (0.3 µm particle)	0.02%
Surface	Smooth

\* Herzberg test at 10 cm hydrostatic head



## 934-AH generalized product protocol

Filter a measured quantity of water (usually 100 ml) through a circle of 934-AH that has been previously washed, dried, and weighed (Fig 3).

### 934-AH RTU Glass Microfiber Filters

#### The filter you trust most just got better.

The 934-AH RTU meets the requirements under the latest edition of Standard Methods 2540D for determining total suspended solids in water. It comes pre-washed and pre-weighed. This eliminates filter pre-treating and streamlines your laboratory process. Each RTU filter is packaged within an aluminum pan with the filter weight clearly printed on a heat resistant label (see Fig. 4).

#### 1. Determination of suspended solids

After filtering, wash solids retained on the 934-AH circle with distilled water and dry at specified temperature (e.g., 103°C–105°C) for a specified time (30 min to 1 h).

#### 2. Determination of total dissolved solids

Total dissolved solids (filterable residue) is that part which passes through a standard glass microfiber filter. The filtrate is collected, evaporated, and dried to a constant weight at a specified temperature and time (drying temperature influences results because weight losses due to volatilization or decomposition and gains due to oxidation all depend on temperature and period of heating).

#### 3. Determination of total solids

Total solids is sum of (1) nonfilterable and (2) filterable solids.

#### 4. Determination of total volatile solids

Total volatile solids is a measure of how much organic matter is present in effluent and an indication of bacterial levels. Measurement can be taken immediately after the total suspended solids determination. Residue on the 934-AH is ignited in a muffle furnace at 550°C. The loss in weight on ignition is measured in mg/l and indicates total volatile solids.



**Fig 3.** Whatman 934-AH RTU (Ready to use) Glass Microfiber filters for determining total suspended solids in water.

### Ordering information

Catalog number	Dimensions (mm)
1827-021	21
1827-024	24
1827-025	25
1827-028	28
1827-030	30
1827-032	32
1827-035	35
1827-037	37
1827-042	42.5
1827-047	47
1827-055	55
1827-070	70
1827-090	90
1827-105	105
1827-110	110
1827-125	125
1827-150	150
1827-185	185
1827-240	240
1827-320	320
9907-042 (RTU format)	42.5
9907-047 (RTU format)	47
9907-055 (RTU format)	55
9907-090 (RTU format)	90
9907-110 (RTU format)	110

GE Healthcare Bio-Sciences Corp.  
800 Centennial Avenue  
Piscataway, NJ 08855-1327 USA

For more information, visit  
[www.gelifesciences.com/whatman](http://www.gelifesciences.com/whatman)



GE, imagination at work and GE monogram are trademarks of General Electric Company. Whatman and 934-AH are trademarks of GE Healthcare companies. All third party trademarks are the property of their respective owners. © 2001-2009 General Electric Company – All rights reserved. Previously published Jan. 2001.

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare Representative for the most current information.