


# Whatman Provides Complete Solutions

Speed your microbiological research with our full range of innovative separations products.

Whatman®  
Part of GE Healthcare  
Visit [www.whatman.com](http://www.whatman.com) or  
call 1-800-WHATMAN for more  
information or to place an order.



## Membranes

## GD/X®

## Mini-UniPrep™

## Puradisc™

## Monitors & Media

### Advantage:

Extensive range of sizes and micron ratings

Easier filtration of high particulate solutions

Prepare HPLC samples three times faster

Efficient filtration of samples up to 100 mL

Pre-sterilized to simplify the testing process

### Feature:

- Consistent high quality and superior performance
- Media include cellulose, nylon, polypropylene and PTFE
- Extreme precision in pore density and pore size distribution
- Suitable for air filtration, aqueous liquids and solvents

- A high quality graded density glass microfiber prefilter
- Available for 13 mm and 25 mm diameter syringe filters
- Sterile and non-sterile options
- Increases volume of sample filtered compared with conventional filters

- Syringeless all-in-one filtration process
- Wide range of membranes to choose from
- Slit septa cap available for automation
- Amber vials available to protect against UV damage

- Available in 4 mm, 13 mm and 25 mm
- Tube tip outlet for 4 mm and 13 mm allows sample to be dispensed accurately in a micro-vial, eliminating air lock
- 4 mm ideal for small volume sample preparation
- 13 mm and 25 mm ideal for HPLC sample preparation

- Ready-to-use filtering units save 70% in time
- Fast and safe—no flaming required
- Unit easily converts to Petri dish for culturing
- Wide range of media for culturing a variety of samples