

Media for a Complete Range of Diagnostic Applications

No matter which diagnostic applications you target, Whatman has the perfect solid-phase media to help expedite your manufacturing and development efforts. Whatman components, in compliance with quality management systems, guarantee product performance, consistency and reliability.

The solutions in this brochure are organized by application. To determine the ones that are right for you, use the reference chart below.

WHATMAN DIAGNOSTIC MEDIA SELECTION CHART

Application Grade	Sample Collection	Sample Transport	Sample Preparation	Sample Wicking	Blood Separation	Conjugate Release	Lateral Flow Membrane	Absorbents	Flow Through Test	Dipsticks
CF1				✓						
CF2	✓									
CF3	✓			✓				✓		✓
CF4				✓					✓	✓
CF5	✓							✓	✓	
CF6				✓				✓	✓	
CF7	✓							✓	✓	
BFC180	✓	✓								
VF1			✓		✓				✓	
VF2			✓		✓				✓	
MF1			✓		✓					
LF1			✓		✓					
GF/B			✓	✓					✓	
GF/F			✓	✓						
GF/AVA			✓	✓					✓	
GF/DVA			✓	✓					✓	
Standard 14				✓		✓				
Standard 17				✓		✓				
Rapid 24				✓		✓				
Rapid 27				✓		✓				
PuraBind A							✓		✓	
PuraBind G							✓		✓	
Immunopore							✓			
Anopore			✓						✓	
Track Etched			✓		✓				✓	
FTA	✓	✓	✓							