

## Typical Properties - Qualitative Wet Strengthened Filter Papers

Grade	Properties	Thickness (mm)	Filtration time to Herzberg (s)	Weight (g/m <sup>2</sup> )	Retention- range** (µm)
2294	fast, wet strength, thick	1.5	55	550	8 – 15
2589 a	medium fast, wet strength	0.45	120	200	6 – 12
858	medium fast, wet strength	0.21	110	75	7-12
860	medium fast, wet strength	0.18	120	75	7-12
0903	medium to slow, wet strength	0.15	350	65	4-7
0905	very fast	0.24	40	65	12-25
Shark Skin	medium to slow, wet strength, thin	0.17	155	44	8 – 12

\*\* Approximate values